

affigene® product list

affigene® trender - quantitation with real-time PCR detection

Product	Code No.	No. of tests	PCR platform (CE, IVD)	Sample preparation (CE, IVD)
affigene® BKV trender	71001084	48	Mx, RG, iC, ABI	affigene®, eM
affigene® CMV trender	71000940	48	Mx, RG, iC, ABI	affigene®, eM
affigene® EBV trender	71000938	48	Mx, RG*, iC, ABI	affigene®, eM
affigene® HBV trender	71000939	48	Mx, RG, iC, ABI	affigene®, eM

affigene® tracer - qualitative determinations with real-time PCR detection

affigene® Aspergillus tracer	71000935	48	Mx, RG*, iC, ABI	
affigene® Cp/Mp tracer	71001306	48	Mx, RG, iC, ABI	affigene®, eM
affigene® HSV 1/2 tracer	71001063	48	Mx, RG, iC, ABI	affigene®, eM
affigene® VZV tracer	71001306	48	Mx, RG, iC, ABI	affigene®, eM

Sample preparation for affigene® trender/tracer products

affigene® DNA extraction	71000942
--------------------------	----------

RG: Rotor-Gene® 3000 and 6000 (RG*: Rotor-Gene® 3000) (Qiagen/Corbett Research)
Mx: Mx3000P® and Mx3005P® (Agilent/Stratagene)
iC: iCycler™ IQ/IQ5 (BioRad)
ABI: ABI 7300 (Applied Biosystems®)
affigene® : affigene® DNA extraction (Cepheid AB)
eM: NucliSENS® easyMAG® (bioMérieux)

